

Thank you for purchasing this Croydex product. Please follow the fixing and care instructions to ensure your grab bar is fitted safely and it retains its high quality finish.

Fittings instruction for the following products:

AP5010** - 300mm Stainless Steel Grab Bar with Concealed Fixings
AP5011** - 450mm Stainless Steel Grab Bar with Concealed Fixings
AP5012** - 600mm Stainless Steel Grab Bar with Concealed Fixings
AP5013** - Angled Stainless Steel Grab Bar with Concealed Fixings
AP500541 - 300mm Stainless Steel Grab Bar with Anti-Slip Grip
AP500741 - 600mm Stainless Steel Grab Bar with Anti-Slip Grip

Tools Required

Pencil, Drill, Masonry drill bit appropriate to fixings, Ceramic drill bit (optional) appropriate to fixings, Spirit Level, Screwdriver appropriate to fixings.

Please ensure that suitable fixings are purchased that are compatible with wall construction materials

KEEP DIY TOOLS OUT OF THE REACH OF CHILDREN

Safety and Care Advice

1. It is recommended that two people fit this product to ensure that it is fitted safely.
2. Always take care when using power tools, particularly in the bathroom.
3. Always check for hidden cables and pipework before drilling and take extreme care if there is any water in the working area.
4. It is advisable to use a residual current circuit breaker (RCCB) when using mains powered tools.
5. Always wear suitable eye protection when drilling.
6. If the product is to be fitted on a ceramic tiled wall, a ceramic drill bit should be used.
7. Always ensure that the drill hole passes through the central tiled area rather than through the grouted area. To prevent damage to the tile, mask the area around the hole with tape before drilling.
8. Do not strike the product with hard or sharp objects.

Product Location

Grab and support bars are only as strong as the wall to which they are fixed and the fixings that are used.

For different wall types you will need to ensure that you are using the correct type of fixing for the material of the wall.

The addition of backboards on the outside of the wall or appropriate strengthening of internal structures may be required to provide adequate wall strength.

If the wall is made of lightweight thermal blocks, aerated hollow brickwork or prefabricated partitions even the most appropriate fixings may not be able to withstand the loads which can be suddenly applied to grab bars, hinged support bars and seats.

Fitting Instruction



Fig.1




Fig.2

Fixings are shown for reference only and are not supplied with the product.

1. Separate the wall bracket covers from the wall brackets by twisting the covers in an anti-clockwise direction.
2. Place the grab bar in position on the wall and using the hole centres of each wall bracket as a template, mark the screw fixing positions on the wall surface. If fitting the product in a vertical or horizontal orientation it is important to ensure the grab bar is level. (See Fig.2)
3. Drill the holes in the marked positions using a masonry drill. (See Fig.2) If the product is to be fitted onto a ceramic tiled wall, a ceramic drill bit should be used.
4. Insert wall plugs flush to the surface of the wall. If fitting to a tiled surface, the wall plugs should be inserted below the surface of the tile to avoid cracking.
5. Insert the screws through the holes of the wall brackets and screw to the wall. (See Fig 2)
6. Slide the wall bracket covers back onto the wall brackets and twist in a clockwise direction to lock the covers in position. (See Fig.1)

Additional Information

- **This product must only be fitted to solid brick, concrete block walls or similarly strong materials offering appropriate support for the maximum safe working load.**
- Maximum safe working load for the Grab Bar is 100kg.

Product Care Advice

- To retain the best quality finish, clean product regularly with a soft damp cloth.
- Do not use abrasive or chemical cleaners, as these will damage the product.

Please retain this leaflet for future reference.

Croydex, Central Way, Andover, Hants, SP10 5AW, UK
Tel: +44 (0) 1264 365881. Fax: +44 (0) 1264 356437
Website: www.croydex.com E-mail: info@croydex.co.uk.